smart Löffelmaterial

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Instructions for use

Patient target group
All persons for whom individual trays are required for treatment.

Intended user
Dentist

Composition Qualitative data:
Oligomeric urethane acrylate, vinyl ester resin, 2,2 dimethoxy-1,2-diphenylethane-1-one or with 2,2 dimethoxy-2-phenylacetophenone, (+)-2,3-dioxo-1,7,7-trimethylbicyclo(2,2,1)heptan-2,3-dione, Diethyl-2,4,6-trimethylbenzoylphosphine oxide, 2-hydroxy-2-methyl-1-phenylpropane-1-one, 2-ethylhexyl-4-(dimethylamino)benzoate, 2,2-diethoxyacetophenone, soda-line glass, amorphous silica, dyes (Pigment Red 146 = C12485 / Pigment Blue 29 = C177007 / D&C Red 30 = C173360 / Pigment Blue C174160 / Pigment Violet 19 = C173900), flavorings (Optamint).

Quantitative data:
In 100 g there are 19 g acrylates/methacrylates and 80 g glass.

Intended purpose of the product
The product is a light-curing tray material in the form of prefabricated templates for the fabrication of torsion-resistant customised impression trays. Alternatively available thinner material (base plates) can be used for the fabrication of bite templates, bite registrations and trays.

Indications
Custom impression tray

Contraindications
The material contains methacrylates. Do not use in patients with an allergy to methacrylates or other substances contained in the product.

Performance characteristics of the product
Smart Löffelmaterial polymerises with light to form a firm and compact polymer.

Specific precautions
Avoid any contact with eyes and skin and the uncured material, i.e. also with the innermost layer. Due to the dust development during the brushing, use a vacuum extraction.

Preparation/Manufacture
Mark the outline of the future tray on the model.

Block out the undercuts on the model.

Use a spacer for the impression material, block out the entire model (e.g. with a sheet of wax) and cut out steps.

If the plaster is very dry, soak the model briefly in water or apply dry isolation or an alginate isolation before use.

Remove the baseplate from the box and close it.

Once the model has been blocked out, adapt the baseplate to it and cut to shape using, e.g. a scalpel.

Use the remaining material to form a handle and attach it to the tray.

Place the model and tray in a light curing unit and cure without vacuum.

Remove the tray from the model and polymerize once again from the underside.

If the tray has been cured without a vacuum, use alcohol to remove the smear layer.

Times:
Working at daylight 20 minutes

Polymerisation in an ultraviolet or LED light curing unit, depending on the unit in use: The time of the tray, on the model 3-5 minutes

Trimming:
Trim the custom tray with tungsten carbide burs and smooth it with sandpaper.

Hygienic
Immersion impressions taken with smart Löffelmaterial custom trays in a disinfectant solution based on glutaraldehyde.

The finished medical product is used by the dentist, e.g. for individual trays with impression material.

Caution
The product has been developed for use in the dental sector and must be used according to the directions. The manufacturer accepts no liability for damage that results from use other than the above.

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Side effects
Adverse side effects are not known.

Interactions
Interaction of the product with other medical devices, medicinal products or other substances is not known.

Storage
The Product should be stored between 0°C - 22°C.

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