

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878 Issue date: 11/29/2023 Revision date: 11/29/2023 Supersedes version of: 11/28/2022 Version: 6.00

Email competent person

sds@kft.de

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier	
Trade name Article number	 smartwipes Wet / smartwipes wet&easy REF 125220 smartwipes Wet neutral; REF 203458 smartwipes Wet lemon; REF 203459 smartwipes Wet ocean; REF 203460 smartwipes Wet tropic; REF 203461 smartwipes Wet bubble gum; REF 214375 smartwipes wet&easy
1.2. Relevant identified uses of the	substance or mixture and uses advised against
1.2.1. Relevant identified uses	
Use of the substance/mixture	: Medical device Germicide Cleaning agent
1.2.2. Lloss advised against	

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Alfred Becht GmbH Carl-Zeiss-Str. 16 P.O. Box 1145 77656 Offenburg T +49 781 60586-0 - F +49 781 60586-40

1.4. Emergency telephone number

Emergency number

: Poison Information Center Freiburg + 49 761 19240

SECTION 2: Hazards identification

Classification according to Regulation (EC) No. 1272/2008 [CLP]		
Flammable liquids, Category 2	H225	
Serious eye damage/eye irritation, Category 1		
Full text of H- and EUH-statements: see section 16		

Adverse physicochemical, human health and environmental effects

Highly flammable liquid and vapour. Causes serious eye damage.

2.2. Label elements

Labelling according to Regulation (EC) No.	1272/2008	[CLP]
Hazard pictograms (CLP)	:	$\mathbf{\Lambda}$

GHS02 GHS05
: Danger
: n-propanol
: H225 - Highly flammable liquid and vapour.
H318 - Causes serious eye damage.
: P210 - Keep away from heat, hot surfaces, sparks, open flames. No smoking.
P280 - Wear eye protection.
P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

contact lenses, if present and easy to do. Continue rinsing. P310 - Immediately call a POISON CENTER, a doctor.

2.3. Other hazards

PBT: not relevant – no registration required

vPvB: not relevant - no registration required

Contains no PBT and/or vPvB substances ≥ 0.1% assessed in accordance with REACH Annex XIII

Component	
n-propanol (71-23-8)	This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII
butanone (78-93-3)	This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII
Ethanol (64-17-5)	This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII
didecyldimethylammonium chloride (7173-51-5)	This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures Comments

: Cleaning agent

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Ethanol substance with national workplace exposure limit(s) (DE)	CAS-No.: 64-17-5 EC-No.: 200-578-6 EC Index-No.: 603-002-00-5	≥ 25 – < 50	Flam. Liq. 2, H225 Eye Irrit. 2, H319
n-propanol	CAS-No.: 71-23-8 EC-No.: 200-746-9 EC Index-No.: 603-003-00-0	≥5-<10	Flam. Liq. 2, H225 Eye Dam. 1, H318 STOT SE 3, H336
butanone substance with national workplace exposure limit(s) (DE); substance with a Community workplace exposure limit	CAS-No.: 78-93-3 EC-No.: 201-159-0 EC Index-No.: 606-002-00-3	≥ 0.25 – < 1	Flam. Liq. 2, H225 Eye Irrit. 2, H319 STOT SE 3, H336 EUH066
didecyldimethylammonium chloride	CAS-No.: 7173-51-5 EC-No.: 230-525-2 EC Index-No.: 612-131-00-6	< 0.1	Acute Tox. 3 (Oral), H301 (ATE=264 mg/kg bodyweight) Skin Corr. 1B, H314 Eye Dam. 1, H318 Aquatic Acute 1, H400 (M=10) Aquatic Chronic 2, H411 EUH071

Specific concentration limits:		
Name	Product identifier	Specific concentration limits (%)
Ethanol	CAS-No.: 64-17-5 EC-No.: 200-578-6 EC Index-No.: 603-002-00-5	(50 ≤ C < 100) Eye Irrit. 2, H319

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

Full text of H- and EUH-statements: see section 16

SECTION 4: First aid measures		
4.1. Description of first aid measures		
First-aid measures general	: When in doubt or if symptoms are observed, get medical advice.	
First-aid measures after inhalation	: Remove person to fresh air and keep comfortable for breathing.	
First-aid measures after skin contact	: Rinse skin with water/shower. Take off immediately all contaminated clothing.	
First-aid measures after eye contact	: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Call a physician immediately.	
First-aid measures after ingestion	: Call a poison center or a doctor if you feel unwell.	
4.2. Most important symptoms and effects, I	both acute and delayed	
Symptoms/effects after eye contact	: Serious damage to eyes.	

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures		
5.1. Extinguishing media		
Suitable extinguishing media Unsuitable extinguishing media	 Use extinguishing media appropriate for surrounding fire. Water spray. Dry powder. Foam. Carbon dioxide. Strong water jet. 	
5.2. Special hazards arising from the substance or mixture		
Fire hazard Explosion hazard Hazardous decomposition products in case of fire	 Highly flammable liquid and vapour. Explosive vapour/air mixtures may be formed. Toxic fumes may be released. Carbon dioxide (CO2). Carbon monoxide. 	
5.3. Advice for firefighters		
Firefighting instructions Protection during firefighting Other information	 Protect container with water spray. Do not attempt to take action without suitable protective equipment. Self-contained breathing apparatus. Complete protective clothing. Do not allow run-off from fire fighting to enter drains or water courses. Cool closed containers exposed to fire with water spray. Disposal must be done according to official 	

SECTION 6: Accidental release measures	
6.1. Personal precautions, protective equ	uipment and emergency procedures
General measures	: Remove all sources of ignition. Avoid contact with skin and eyes. Ensure adequate air ventilation. Do not breathe gas/vapour/aerosol.
6.1.1. For non-emergency personnel	
Emergency procedures	: Ventilate spillage area. No open flames, no sparks, and no smoking. Avoid contact with skin and eyes.
6.1.2. For emergency responders	
Protective equipment	: Do not attempt to take action without suitable protective equipment. For further information refer to section 8: "Exposure controls/personal protection".
6.2. Environmental precautions	

Avoid release to the environment. Prevent entry to sewers and public waters. Avoid sub-soil penetration.

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

6.3. Methods and material for containment and cleaning up		
Methods for cleaning up	: Cover spill with non combustible material, e.g.: sand/earth. Take up mechanically (sweeping, shovelling) and collect in suitable container for disposal.	
Other information	: Dispose of materials or solid residues at an authorized site. Disposal must be done according to official regulations.	

6.4. Reference to other sections

Information for safe handling. See section 7. Concerning personal protective equipment to use, see section 8. For further information refer to section 13.

SECTION 7: Handling and storage			
7.1. Precautions for safe handling			
Additional hazards when processed Precautions for safe handling Hygiene measures	 In use, may form flammable vapour-air mixture. Ensure good ventilation of the work station. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Take precautionary measures against static discharge. Flammable vapours may accumulate in the container. Use only non-sparking tools. Use explosion-proof equipment. Wear personal protective equipment. Avoid contact with skin and eyes. Immediately remove contaminated or damp clothing. Do not eat, drink or smoke when using this product. Always wash hands after handling the product. 		
7.2. Conditions for safe storage, including any incompatibilities			
Technical measures	: Ground/bond container and receiving equipment.		
Storage conditions	: Protect against frost. Store in a well-ventilated place. Keep cool. Keep container tightly closed.		
Heat and ignition sources	: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep away from heat and direct sunlight.		
Information about storage in one common storage facility	: Keep away from food, drink and animal feeding stuffs. Store away from Acids, oxidizing materials.		

7.3. Specific end use(s)

Follow the directions!.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

8.1.1 National occupational exposure and biological limit values

butanone (78-93-3)		
EU - Indicative Occupational Exposure Limit (IOEL)		
Local name	Butanone	
IOEL TWA	600 mg/m³	
IOEL TWA [ppm]	200 ppm	
IOEL STEL	900 mg/m³	
IOEL STEL [ppm]	300 ppm	
Regulatory reference	COMMISSION DIRECTIVE 2000/39/EC	
Germany - Occupational Exposure Limits (TRGS 900)		
Local name	Butanon	
AGW (OEL TWA) [1]	600 mg/m³	
AGW (OEL TWA) [2]	200 ppm	

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

butanone (78-93-3)		
Peak exposure limitation factor	1(I)	
Remark	DFG - Senatskommission zur Prüfung gesundheitsschädlicher Arbeitsstoffe der DFG (MAK-Kommission); EU - Europäische Union (Von der EU wurde ein Luftgrenzwert festgelegt: Abweichungen bei Wert und Spitzenbegrenzung sind möglich); H - hautresorptiv; Y - Ein Risiko der Fruchtschädigung braucht bei Einhaltung des Arbeitsplatzgrenzwertes und des biologischen Grenzwertes (BGW) nicht befürchtet zu werden	
Regulatory reference	TRGS900	
Germany - Biological limit values (TRGS 903)		
Local name	2-Butanon (Methylethylketon)	
Biological limit value	2 mg/l Parameter: 2-Butanon - Untersuchungsmaterial: U = Urin - Probenahmezeitpunkt: b) Expositionsende, bzw. Schichtende - Festlegung/Begründung: 05/2015 DFG	
Remark	U, b, 05/2015, DFG	
Regulatory reference	TRGS 903	
Ethanol (64-17-5)		
Germany - Occupational Exposure Limits (TRGS 900)		
Local name	Ethanol	
AGW (OEL TWA) [1]	380 mg/m ³	
AGW (OEL TWA) [2]	200 ppm	
Peak exposure limitation factor	4(II)	
Remark	DFG - Senatskommission zur Prüfung gesundheitsschädlicher Arbeitsstoffe der DFG (MAK-Kommission); Y - Ein Risiko der Fruchtschädigung braucht bei Einhaltung des Arbeitsplatzgrenzwertes und des biologischen Grenzwertes (BGW) nicht befürchtet zu werden	
Regulatory reference	TRGS900	

8.1.2. Recommended monitoring procedures

No additional information available

8.1.3. Air contaminants formed

No additional information available

8.1.4. DNEL and PNEC

n-propanol (71-23-8)		
DNEL/DMEL (Workers)		
Acute - systemic effects, inhalation	1037 mg/m ³	
DNEL/DMEL (General population)		
Acute - systemic effects, inhalation	518 mg/m³	
butanone (78-93-3)		
DNEL/DMEL (Workers)		
Acute - systemic effects, inhalation	900 mg/m³	
Long-term - systemic effects, dermal	1161 mg/kg bodyweight/day	
Long-term - systemic effects, inhalation	600 mg/m³	
DNEL/DMEL (General population)		
Acute - systemic effects, inhalation	450 mg/m ³	

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

butanone (78-93-3)		
Long-term - systemic effects,oral	31 mg/kg bodyweight/day	
Long-term - systemic effects, inhalation	106 mg/m ³	
Long-term - systemic effects, dermal	412 mg/kg bodyweight/day	
Ethanol (64-17-5)		
DNEL/DMEL (Workers)		
Long-term - systemic effects, dermal	8238 mg/kg bodyweight/day	
Long-term - systemic effects, inhalation	380 mg/m³	
DNEL/DMEL (General population)		
Long-term - systemic effects, inhalation	114 mg/m ³	
PNEC (Water)		
PNEC aqua (freshwater)	0.96 mg/l	
PNEC aqua (marine water)	0.79 mg/l	
PNEC aqua (intermittent, freshwater)	2.75 mg/l	
PNEC (Sediment)		
PNEC sediment (freshwater)	3.6 mg/kg dwt	
PNEC sediment (marine water)	2.9 mg/kg dwt	
PNEC (Soil)		
PNEC soil	0.63 mg/kg dwt	
PNEC (Oral)		
PNEC oral (secondary poisoning)	380 mg/kg food	
PNEC (STP)		
PNEC sewage treatment plant	580 mg/l	
didecyldimethylammonium chloride (7173-51-5)		
PNEC (Water)		
PNEC aqua (freshwater)	0.0011 mg/l	
PNEC aqua (marine water)	0.00011 mg/l	
PNEC aqua (intermittent, freshwater)	0.00021 mg/l	
PNEC aqua (intermittent, marine water)	0.000021 mg/l	
PNEC (Sediment)		
PNEC sediment (freshwater)	61.86 mg/kg dwt	
PNEC sediment (marine water)	6.186 mg/kg dwt	
PNEC (Soil)		
PNEC soil	1.4 mg/kg dwt	
PNEC (STP)		
PNEC sewage treatment plant	0.14 mg/l	

8.1.5. Control banding

No additional information available

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

8.2. Exposure controls

8.2.1. Appropriate engineering controls

Appropriate engineering controls:

Ensure good ventilation of the work station. Emergency eye wash fountains should be available in the immediate vicinity of any potential exposure.

8.2.2. Personal protection equipment

8.2.2.1. Eye and face protection

Eye protection:

Wear closed safety glasses. ISO 16321-1. Emergency eye wash fountains should be available in the immediate vicinity of any potential exposure

8.2.2.2. Skin protection

Skin and body protection:

Wear suitable protective clothing. EN 13034. EN ISO 13688

Hand protection:

In case of repeated or prolonged contact wear gloves. Chemically resistant protective gloves. Butyl rubber. Neoprene. Polyethylene. ISO 374-1. Choosing the proper glove is a decision that depends not only on the type of material, but also on other quality features, which differ for each manufacturer. Please follow the instructions related to the permeability and the penetration time provided by the manufacturer. Gloves must be replaced after each use and whenever signs of wear or perforation appear

8.2.2.3. Respiratory protection

Respiratory protection:

In case of insufficient ventilation, wear suitable respiratory equipment. EN 143. Short term exposure. Breathing apparatus with filter. A-P2. Breathing equipment is only to be used in order to handle the residual risk of short term jobs if all other risk minimizing measures have been carried out e.g. retention and/or local exhaust. For details on conditions of use and maximum use concentrations, see DGUV Regulation 112-190 - Use of respiratory protective equipment.

8.2.2.4. Thermal hazards

No additional information available

8.2.3. Environmental exposure controls

Environmental exposure controls:

Avoid release to the environment.

Other information:

Do not eat, drink or smoke when using this product. Avoid contact with skin and eyes. Always wash hands after handling the product. Apply emollient cream.

SECTION 9: Physical and chemical properties		
9.1. Information on basic physical and chemical properties		
Physical state	: Solid	
Colour Appearance	: white. : Inert medium impregnated with: Germicide.	
Odour	: alcohol odour.	
Odour threshold	: Not available	
Melting point	: Not applicable	
Freezing point	: Not available	
Boiling point	: Not available	
Flammability	: Not applicable	
Explosive properties	: Product is not explosive. Explosive vapour/air mixtures may be formed.	
Oxidising properties	: Non oxidizing.	
Lower explosion limit	: Not applicable	
Upper explosion limit	: Not applicable	
Flash point	: Not applicable	
Auto-ignition temperature	: Not applicable	
Decomposition temperature	: Not available	
рН	: <10 (Solution)	
pH solution	: Not available	
Viscosity, kinematic	: Not applicable	

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

Solubility Partition coefficient n-octanol/water (Log Kow) Partition coefficient n-octanol/water (Log Pow) Vapour pressure Vapour pressure at 50°C Density Relative density Relative vapour density at 20°C	 Water: Miscible (Solution) Not available Not applicable Not available Not available 0.91 g/cm³ (Solution) Not available Not available Not available Not available Not available
Relative vapour density at 20°C Particle size	: Not applicable : Not available

9.2. Other information

9.2.1. Information with regard to physical hazard classes

No additional information available

9.2.2. Other safety characteristics

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

Highly flammable liquid and vapour.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use. Explosive vapour/air mixtures may be formed.

10.4. Conditions to avoid

Avoid contact with hot surfaces. Heat. No flames, no sparks. Eliminate all sources of ignition.

10.5. Incompatible materials

No additional information available

10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced.

11.1. Information on hazard classes	s as defined in Regulation (EC) No 1272/2008
Acute toxicity (oral) Acute toxicity (dermal) Acute toxicity (inhalation)	 Not classified (Based on available data, the classification criteria are not met) Not classified (Based on available data, the classification criteria are not met) Not classified (Based on available data, the classification criteria are not met)
didecyldimethylammonium chloride (7173-51-5)	
LD50 oral rat	264 mg/kg bodyweight (female; (OECD 401 method))
LD50 dermal rat	> 2000 mg/kg bodyweight (OECD 402 method)
Skin corrosion/irritation	: Not classified (Based on available data, the classification criteria are not met) pH: < 10 (Solution)
Serious eye damage/irritation	: Causes serious eye damage. pH: < 10 (Solution)
Respiratory or skin sensitisation	: Not classified (Based on available data, the classification criteria are not met)
Germ cell mutagenicity	: Not classified (Based on available data, the classification criteria are not met)
Carcinogenicity	: Not classified (Based on available data, the classification criteria are not met)
Reproductive toxicity	: Not classified (Based on available data, the classification criteria are not met)

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

STOT-single exposure :	Not classified (Based on available data, the classification criteria are not met)
n-propanol (71-23-8)	
STOT-single exposure	May cause drowsiness or dizziness.
butanone (78-93-3)	
STOT-single exposure	May cause drowsiness or dizziness.
1 1	Not classified (Based on available data, the classification criteria are not met) Not classified (Not relevant)

11.2. Information on other hazards

No additional information available

SECTION 12: Ecological information		
12.1. Toxicity		
	Not classified (Based on available data, the classification criteria are not met)	
(acute) Hazardous to the aquatic environment, long-term : (chronic)	Not classified (Based on available data, the classification criteria are not met)	
didecyldimethylammonium chloride (7173-51-	5)	
LC50 - Fish [1]	0.49 mg/l (96 h; Brachydanio rerio (zebra-fish); (OECD 203 method))	
EC50 - Crustacea [1]	≈ 0.057 mg/l (48 h; Daphnia magna; (OECD 202 method))	
ErC50 algae	0.062 mg/l (72 h; Pseudokirchnerella subcapitata (OECD 201 method))	
NOEC chronic crustacea	0.021 mg/l (21 d; Daphnia magna; (OECD 211 method))	
NOEC chronic algae	0.013 mg/l (OECD 201 method)	
12.2. Persistence and degradability		
smartwipes Wet / smartwipes wet&easy		
Persistence and degradability	The product has not been tested.	
n-propanol (71-23-8)		
Persistence and degradability	Readily biodegradable.	
Biodegradation	75 % (20 d)	
butanone (78-93-3)		
Persistence and degradability	Readily biodegradable.	
Biodegradation	98 % (28 d; (OECD 301D method))	
Ethanol		
Biodegradation	84 % (20 d; Literature data)	
didecyldimethylammonium chloride (7173-51-5)		
Persistence and degradability	Readily biodegradable.	
Biodegradation	69 % (28d)	
12.3. Bioaccumulative potential		
smartwipes Wet / smartwipes wet&easy		
Partition coefficient n-octanol/water (Log Pow)	Not applicable	

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

smartwipes Wet / smartwipes wet&easy		
Bioaccumulative potential	The product has not been tested.	
n-propanol (71-23-8)		
Partition coefficient n-octanol/water (Log Pow)	0.2 (25 °C; pH 7; (OECD 117 method))	
butanone (78-93-3)		
Partition coefficient n-octanol/water (Log Pow)	0.3 (40 °C; pH 7; (OECD 117 method))	
Bioaccumulative potential	Bioaccumulation unlikely.	
Ethanol		
Partition coefficient n-octanol/water (Log Pow)	-0.35 (24 °C; pH 7,4; (OECD 107 method))	
didecyldimethylammonium chloride (7173-51	-5)	
Partition coefficient n-octanol/water (Log Pow)	2.59 (20 °C; (OECD 105 method))	
12.4. Mobility in soil		
smartwipes Wet / smartwipes wet&easy		
Ecology - soil	The product has not been tested.	
n-propanol (71-23-8)		
Organic Carbon Normalized Adsorption Coefficient (Log Koc)	0.633 (Quantitative structure-activity relationship (QSAR))	
Ethanol		
Surface tension	22.31 mN/m (20 °C)	
Organic Carbon Normalized Adsorption Coefficient (Log Koc)	≈ 0.2 (Handbook and/or Scientific papers)	
didecyldimethylammonium chloride (7173-51-5)		
Surface tension	25.82 mN/m (OECD 115 method)	
12.5. Results of PBT and vPvB assessment		
smartwipes Wet / smartwipes wet&easy		
PBT: not relevant – no registration required		
vPvB: not relevant – no registration required		
Component		
n-propanol (71-23-8)	This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII	
butanone (78-93-3)	This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII	
Ethanol (64-17-5)	This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII	
didecyldimethylammonium chloride (7173-51-5)	This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII	
	•	

12.6. Endocrine disrupting properties

No additional information available

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

12.7. Other adverse effects

No additional information available

SECTION 13: Disposal considerations	
13.1. Waste treatment methods	
Waste treatment methods Sewage disposal recommendations Product/Packaging disposal recommendations Additional information European List of Waste (LoW, EC 2000/532) HP Code	 European waste catalogue. Disposal must be done according to official regulations. Do not allow into drains or water courses. Do not dispose of with domestic waste. Flammable vapours may accumulate in the container. 07 01 01* - aqueous washing liquids and mother liquors HP3 - "Flammable:" flammable liquid waste: liquid waste having a flash point below 60 °C or waste gas oil, diesel and light heating oils having a flash point > 55 °C and ≤ 75 °C; flammable pyrophoric liquid and solid waste: solid or liquid waste which, even in small quantities, is liable to ignite within five minutes after coming into contact with air; flammable solid waste: solid waste which is readily combustible or may cause or contribute to fire through friction; flammable gaseous waste: gaseous waste which is flammable in air at 20 °C and a standard pressure of 101.3 kPa; water reactive waste: waste which, in contact with water, emits flammable gases in dangerous quantities; other flammable waste: flammable aerosols, flammable self-heating waste, flammable organic peroxides and flammable self-reactive waste. HP4 - "Irritant – skin irritation and eye damage:" waste which on application can cause skin irritation or damage to the eye.

SECTION 14: Transport information

ADR	IMDG	ΙΑΤΑ	ADN	RID
14.1. UN number or ID n	umber			
UN 3175	UN 3175	UN 3175	UN 3175	UN 3175
14.2. UN proper shippin	g name			
SOLIDS CONTAINING FLAMMABLE LIQUID, N.O.S. (ethanol ; propan-1- ol)	SOLIDS CONTAINING FLAMMABLE LIQUID, N.O.S. (ethanol ; propan-2- ol)	Solids containing flammable liquid, n.o.s. (ethanol ; propan-2-ol)	SOLIDS CONTAINING FLAMMABLE LIQUID, N.O.S. (ethanol ; propan-2- ol)	SOLIDS CONTAINING FLAMMABLE LIQUID, N.O.S. (ethanol ; propan-2- ol)
Transport document descr	iption			
UN 3175 SOLIDS CONTAINING FLAMMABLE LIQUID, N.O.S. (ethanol ; propan-1- ol), 4.1, II, (E)	UN 3175 SOLIDS CONTAINING FLAMMABLE LIQUID, N.O.S. (ethanol ; propan-2- ol), 4.1, II	UN 3175 Solids containing flammable liquid, n.o.s. (ethanol ; propan-2-ol), 4.1, II	UN 3175 SOLIDS CONTAINING FLAMMABLE LIQUID, N.O.S. (ethanol ; propan-2- ol), 4.1, II	UN 3175 SOLIDS CONTAINING FLAMMABLE LIQUID, N.O.S. (ethanol ; propan-2- ol), 4.1, II
14.3. Transport hazard o	class(es)			
4.1	4.1	4.1	4.1	4.1

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

ADR	IMDG	ΙΑΤΑ	ADN	RID
14.4. Packing group				
II	II	II	II	II
14.5. Environmental haz	ards		1	
Dangerous for the environment: No	Dangerous for the environment: No Marine pollutant: No	Dangerous for the environment: No	Dangerous for the environment: No	Dangerous for the environment: No
No supplementary information	n available			
14.6. Special precautions	s for user			
Overland transport Classification code (ADR) Special provisions (ADR) Limited quantities (ADR) Excepted quantities (ADR) Transport category (ADR) Hazard identification number (Orange plates	: 1kg : E2 : 2	-		
Tunnel restriction code (ADR)	: E			
Transport by seaSpecial provisions (IMDG): 24EmS-No. (Fire): F-EmS-No. (Spillage): S-				
Air transportPCA Excepted quantities (IATA): E2PCA Limited quantities (IATA): Y4PCA limited quantity max net quantity (IATA): 5kPCA packing instructions (IATA): 44PCA max net quantity (IATA): 15CAO max net quantity (IATA): 50Special provisions (IATA): A44		41 9 5 &g &g		
Inland waterway transport Classification code (ADN) Special provisions (ADN) Limited quantities (ADN) Excepted quantities (ADN) Carriage permitted (ADN) Additional requirements/Rema	: 1 k : E2 : B	-	when this substance is carried	t in bulk or without packagin
Rail transport Classification code (RID) Special provisions (RID) Limited quantities (RID) Excepted quantities (RID) Transport category (RID) Hazard identification number (: 1kg : E2 : 2	-		

14.7. Maritime transport in bulk according to IMO instruments

Not applicable

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Other information, restriction and prohibition : Take note of Directive 94/33/EC on the protection of young people at work. regulations

REACH Annex XVII (Restriction List)

EU restriction list (REACH Annex XVII)		
Reference code	ference code Applicable on	
3(a)	smartwipes Wet / smartwipes wet&easy ; butanone ; Ethanol ; n-propanol	
3(b)	smartwipes Wet / smartwipes wet&easy ; butanone ; Ethanol ; n-propanol	
40.	smartwipes Wet / smartwipes wet&easy ; butanone ; Ethanol ; n-propanol	

REACH Annex XIV (Authorisation List)

Contains no substance(s) listed on REACH Annex XIV (Authorisation List)

REACH Candidate List (SVHC)

Contains no substance(s) listed on the REACH Candidate List

PIC Regulation (Prior Informed Consent)

Contains substance(s) listed on the PIC list (Regulation EU 649/2012 concerning the export and import of hazardous chemicals): Didecyldimethylammonium chloride (7173-51-5)

POP Regulation (Persistent Organic Pollutants)

Contains no substance(s) listed on the POP list (Regulation EU 2019/1021 on persistent organic pollutants)

Ozone Regulation (1005/2009)

Contains no substance(s) listed on the Ozone Depletion list (Regulation EU 1005/2009 on substances that deplete the ozone layer)

Seveso Directive (Disaster Risk Reduction)

Seveso III Part I (Categories of dangerous substances)	Qualifying quantity (tonnes)		
	Lower-tier	Upper-tier	
P5c FLAMMABLE LIQUIDS Flammable liquids, Categories 2 or 3 not covered by P5a and P5b	5000	50000	

Explosives Precursors Regulation (2019/1148)

Contains no substance(s) listed on the Explosives Precursors list (Regulation EU 2019/1148 on the marketing and use of explosives precursors)

Drug Precursors Regulation (273/2004)

Contains substance(s) listed on the Drug Precursors list (Regulation EC 273/2004 on the manufacture and the placing on market of certain substances used in the illicit manufacture of narcotic drugs and psychotropic substances)

Name	CN designation	CAS-No.	CN code	Category	Threshold	Annex
Methylethylketone	Butanone	78-93-3	2914 12 00	Category 3		Annex I

15.1.2. National regulations

Germany

Employment restrictions

: Employment prohibitions or restrictions on the protection of young people at work according to § 22 JArbSchG in the case of formation of hazardous substances have to be observed.

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

National Rules and Recommendations	: TRGS 400: Risk Assessment for Activities involving Hazardous Substances.
	TRGS 500: Protective measures.
	TRGS 510: Storage of hazardous substances in non-stationary containers.
	TRGS 520: Construction and operation of collection points and temporary storage for small amounts of hazardous waste.
	TRGS 900: Occupational Exposure Limits.
	TRGS 903: Biological limit values.
Water hazard class (WGK)	: WGK 1, Slightly hazardous to water (Classification according to AwSV, Annex 1).
Storage class (LGK, TRGS 510)	: LGK 4.1B - Flammable solids.
Hazardous Incident Ordinance (12. BImSchV)	: Listed in the 12. BlmSchV (Annex I) under: 1.2.5.3
	- Quantity threshold for operational area under § 1 para. 1
	- Sentence 1 :5000000 kg
	- Sentence 2 :5000000 kg

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

SECTION 16: Other information

Indication of changes:

General revision.

Indication of changes			
Section	Changed item	Change	Comments
	General revision		
15	Storage class (LGK)	Modified	

Abbreviations and acronyms:			
ADN	European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways		
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road		
BCF	Bioconcentration factor		
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008		
EC50	Median effective concentration		
ΙΑΤΑ	International Air Transport Association		
IMDG	International Maritime Dangerous Goods		
LC50	Median lethal concentration		
LD50	Median lethal dose		
РВТ	Persistent Bioaccumulative Toxic		
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006		
RID	Regulations concerning the International Carriage of Dangerous Goods by Rail		
SDS	Safety Data Sheet		
vPvB	Very Persistent and Very Bioaccumulative		
ATE	Acute Toxicity Estimate		
DNEL	Derived-No Effect Level		
DMEL	Derived Minimal Effect level		
IARC	International Agency for Research on Cancer		
LOAEL	Lowest Observed Adverse Effect Level		

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

Abbreviations and acronyms:		
NOAEC	No-Observed Adverse Effect Concentration	
NOAEL	lo-Observed Adverse Effect Level	
NOEC	No-Observed Effect Concentration	
OECD	Organisation for Economic Co-operation and Development	
PNEC	Predicted No-Effect Concentration	
STP	Sewage treatment plant	
TLM	Median Tolerance Limit	

Data sources

: Information provided by the manufacturer. MSDSs of the suppliers. European Chemicals Agency, http://echa.europa.eu/.

Department issuing data specification sheet:

KFT Chemieservice GmbH
Im Leuschnerpark 3
D-64347 Griesheim
Phone: +49 6155-8981-400
Fax: +49 6155 8981-500

SDS Service: +49 6155 8981-522 : Dr. Maximilian Gatterdam

Contact person

Full text of H- and EUH	Full text of H- and EUH-statements:		
Acute Tox. 3 (Oral)	Acute toxicity (oral), Category 3		
Aquatic Acute 1	Hazardous to the aquatic environment – Acute Hazard, Category 1		
Aquatic Chronic 2	Hazardous to the aquatic environment – Chronic Hazard, Category 2		
EUH066	Repeated exposure may cause skin dryness or cracking.		
EUH071	Corrosive to the respiratory tract.		
Eye Dam. 1	Serious eye damage/eye irritation, Category 1		
Eye Irrit. 2	Serious eye damage/eye irritation, Category 2		
Flam. Liq. 2	Flammable liquids, Category 2		
H225	Highly flammable liquid and vapour.		
H301	Toxic if swallowed.		
H314	Causes severe skin burns and eye damage.		
H318	Causes serious eye damage.		
H319	Causes serious eye irritation.		
H336	May cause drowsiness or dizziness.		
H400	Very toxic to aquatic life.		
H411	Toxic to aquatic life with long lasting effects.		
Skin Corr. 1B	Skin corrosion/irritation, Category 1, Sub-Category 1B		
STOT SE 3	Specific target organ toxicity – Single exposure, Category 3, Narcosis		

Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:			
Flam. Liq. 2	H225 Calculation method		
Eye Dam. 1 H318 Calculation method			

KFT SDS EU 01 - Version 23.1

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.